IRSN

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	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 1 of 10

1 Purpose

Acquisition of beam data for the Novalis suitable for Treatment Planning with the BrainSCAN dose algorithm *Pencil Beam*.

2 Scope

This working instruction is valid for all BrainLAB Support Engineers.

3 Definitions and Abbreviations

BPE: Beam Profile Editor
 mMLC: micro Multi Leaf Collimator (Novalis m3)
 NLOUT: Nominal Linac Output
 PBA: Pencil Beam Algorithm
 PDD: Percentage Depth-Dose;
 RFS: Radial Factors
 SFC: Source Function Correction
 SID: Source Isocenter Distance
 SSD: Source Surface Distance
 TMR: Tissue Maximum Ratio

4 Equipment

The following equipment is needed to perform the minimum set of measurements:

- motorized water tank with at least 30x30 cm² base area and a depth of at least 40 cm
- large detector: ionization chamber 0.125 cm³ or 0.3 cm³ (needs to be calibrated for NLOUT measurement)
- medium size detector: pinpoint max. 0.03 cm³
- small detector: either a diode or a diamond detector
- Excel Template provided by BrainLAB

The following equipment is needed for dynamic IMRT installations:

- dynamic mMLC files "M3_1.d01", "M3_5.d01", ... "M3_100.d01"

UNCONTROLLED INFORMATION IF NOT READ ONLINE

Originated:	Approved:	Released:
March 19, 2004	March 19, 2004	April 05, 2004
Váth Albert	Stephan Fröhlich	Thomas Halbwirth
Software Engineer RT	Division Director RT	Technical Support Coordinator

BrainLAB	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 2 of 10

5 Checklist

UNCONTROLLED INFORMATION IF NOT READ ONLINE

Section	Task	Equipment	crucial	Done
6.3.1	Nominal Linac Output (NLOUT)	calibrated ionization chamber	yes	
6.3.2	Leakage for open and closed primary jaws	calibrated ionization chamber	yes	
	Enter NLOUT	Excel template	yes	
	Enter Leakage values	Excel template	yes	
	Enter SSD and Norm. Depth	Excel template	yes	
6.3.3	TMR or PDD Jaw and mMLC field 6x6, 12x12, 18x18, 42x42, 60x60, 80x80, 100x100 [mm ²]	small detector	yes	
	Import TMR or PDD values	Excel template	yes	
6.3.4	Scatter Factors Jaw field 6x6, 12x12, 18x18, 42x42, 60x60, 80x80, 100x100 [mm ²] mMLC field 6x6, 12x12, 18x18, 24x24, 30x30, 36x36, 42x42, 60x60, 80x80, 100x100 [mm ²]	medium size detector	yes	
	Import Scatter Factors	Excel template	yes	
6.3.5	Diagonal Radial Profiles Jaw and mMLC field set to 100x100 mm ² depths: 5, 14, 25, 50, 100, 200, 350,... [mm]	medium size detector	yes	
	Import Diagonal Radial Factors	Excel template	yes	
6.4.1	Dynamic Leaf Shift gap size: 1, 5, 10, 20, 50, 100 [mm]	large detector	dyn. IMRT only	
	Calculate Dynamic Leaf Shift and enter into Beam Profile	Excel Template	dyn. IMRT only	
	Prepare Beam Profile Enter Identifier, Energy, NLOUT, norm. Depth, SSD and Leakage Paste TMR or PDD, Scatter and Radial Factors	Beam Profile Editor	yes	
6.5	Enter default values for SFC and Radiologic Field Correction	Beam Profile Editor	yes	
	Transfer data to BrainLAB support either your local support engineer or send the data to support@brainlab.com		yes	

	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 3 of 10

6 Instructions

6.1 Introduction

The following document describes the measurements that are needed as input data for the *Pencil Beam Algorithm* (PBA).

Minimum measurement requirements for all energies are (Section 6.3):

- Nominal Linac Output (NLOUT) for a field of 100x100 mm²
- Leakage for open and closed primary jaws
- TMR/PDD for 7 jaw settings with mMLC aperture equal to jaw aperture
- Scatter Factors for 7 jaw settings and 10 mMLC square fields, each
- Radial Profiles in diagonal direction for a 100x100 mm² field in 7 depths (to be converted to Radial Factors at BrainLAB HQ)
- Dynamic leaf shift (dynamic IMRT installations, only)

For the measurements a water-phantom of X×Y×Z ≥ 30cm×30cm×40cm is necessary. The beam axis must be vertical, i.e. orthogonal to the water-surface, which should be roughly in 1 m distance from the accelerator source. The highest spatial resolution is needed in the vertical direction for the TMR.

All data has to be entered or copied into the Beam Profile Editor.

Please note, that some additional profile measurements are recommended in the BrainSCAN Software Guide (Appendix D: Pencil Beam Dose Algorithm, Section: Validation). The profiles are intended to be compared with BrainSCAN dose export by the responsible physicist after completion of the beam profile. In order to avoid repeated setup of the water phantom, the measurements may be performed together with the beam profile data.

6.2 Important notes

- Make sure, that the jaw aperture never exceeds 100×100 mm²!
- The accuracy of the BrainSCAN Pencil Beam dose algorithm is directly dependent on the accuracy and the range of the beam data measurement. It has to be assured that the beam data measurement covers the range of field sizes that will be used in later treatment planning. This affects the measurement of the scatter factors, the radial profiles and the PDDs/TMRs.
- For IMRT fields, the relevant field size for the dose algorithm is the one of the corresponding conformal beam and not the field sizes of the subsegments of a beam.
- The number of measurements specified below is completely sufficient for the PBA. If you really want to improve accuracy, perform the measurements with extreme care, repeat them, select the best ones (i.e. lowest noise) and average them. Finer steps in field-size, depth- or radial direction, although not prohibited, will not increase dose accuracy.

UNCONTROLLED INFORMATION IF NOT READ ONLINE

	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 4 of 10

6.3 Measurements for all installations

6.3.1 Measurement of Nominal Linac Output

Use a calibrated detector. This is the only measurement giving an absolute dose value. All the other measurements are relative ones.

- Setup the water phantom, with the isocenter at water surface level (SSD = SID = 1000 mm). Adjust the center of the active detector volume to coincide with the isocenter (level of water surface) and mark this as depth zero (Figure 1)!
- Move the detector to a depth of 50 mm.
- Set the mMLC square field size and the jaw aperture to 100x100 mm².
- Deliver at least 300 MU. The result has to be in Gray (Gy) (electrometer readings, have to be converted).
- Enter the measured value into the field Properties->General Settings->Pencil Beam->"Nominal LINAC Output" together with the "Normalization Depth" of 50 mm and the "Source Surface Distance" of 100 cm.

Remark: It is also possible to use other settings than SSD = 1000 mm and the normalization depth of 50 mm. Especially for higher energies, a larger depth (i. e. 100 mm) may be useful. But all the measurements for NLQUT, leakage, PDD, radial factors in PDD approach and scatter factors have to be done with the same SSD. All the measurements for NLQUT, leakage and scatter factors have to be done at the same depth.


Figure 1: Setup of water phantom and detector for the measurement of the Nominal LINAC Output. The water surface is adjusted at isocenter depth (SSD = SID) and the center of the active detector volume is set at the isocenter (Setup position, depth = 0 mm). For measurement, the detector is moved vertically down to a depth of 50 mm (Measurement).

UNCONTROLLED INFORMATION IF NOT READ ONLINE

	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 5 of 10

UNCONTROLLED INFORMATION IF NOT READ ONLINE

6.3.2 Background leakage for open and closed jaws

For the background leakage the setup is identical to the NLOUT measurement. If the same calibrated detector as for NLOUT is used (recommended!) the NLOUT can be used as reference value. Otherwise the NLOUT measurement has to be repeated with the new sensor to obtain the reference value for the leakage. If the LINAC does not allow an asymmetric setting of the closed jaws, the detector can be moved horizontally in x- and y-directions by 20 mm. In this case the reference value has to be measured at the new position using a square field jaw- and mMLC-aperture of 100 x 100 mm².

- Leave the primary jaws open at 100x100 mm².
- Close the mMLC leaves (asymmetric; the intra-leaf gap should be 50 mm off the isocenter).
- Deliver at least 300 MU.
- Close the primary jaws, too (also asymmetric). Repeat the measurement.
- Enter the dose fractions (relative to NLOUT) in the fields Properties- > Pencil Beam Related Properties- > Multileaf Background Leakage- > "Leakage for Open Jaws" and "Leakage for Closed Jaws".

6.3.3 Measurement of PDD/TMR

Either PDD (percentage depth dose), where the water phantom is fixed and the detector moves along the beam axis or TMR (tissue maximum ratio), where detector is fixed at isocenter and the water surface moves vertically, can be measured. The kind of data entered in the table has to be set using the "Values are ..." - "TMR" or "PDD" buttons. See figure 2 for illustration.

Use a detector with a small active volume to get high spatial resolution. If a diamond detector is used, the orientation of the diamond disk should be horizontal to get the highest resolution in vertical direction. Recalibrate x/y/z-coordinates of the water phantom, when changing the detector. The following steps describe the measurement of PDD values.

- Setup the water phantom in the same way as for the NLOUT measurement (isocenter at water surface level: SSD = SID = 1000 mm).
- Move the jaws as well as the mMLC leaves to form square fields. mMLC and jaw aperture must be equal! The following standard field sizes have to be used for TMR/PDD measurements: 6x6, 12x12, 18x18, 42x42, 60x60, 80x80, 100x100 [mm²].
- Use the water phantom software to measure in depths from 0 to the desired depth. Use a step size of 1mm for 0 - 50 mm depth and a step size of 5 mm for higher depths.
- It is recommended to move the detector from the bottom upward to avoid errors due to the bending of the water surface (capillary effects).
- Take care to measure up to a depth that is equivalent or greater than the maximum depth needed inside the patients body. For cranial treatments 250 mm might be enough, for extracranial cases a range of up to 500 mm or more has to be covered.
- Due to the clearance of your LINAC it will not be possible to measure TMR up to a depth of 300 or 350 mm. If such depths are necessary you have to measure PDD.
- TMR/PDD values can be normalized arbitrarily.


Figure 2: Setup of water phantom and detector for measuring PDD and TMR. In the PDD approach (left sketch), the SSD is fixed and the detector moves vertically along the central beam axis for PDD measurement. Radial factors are measured by moving the detector horizontally in each depth in a diagonal direction with respect to the jaw field. In the TMR approach (right sketch), the detector is fixed at isocenter depth and the TMRs are measured by varying the water level. The RFS are measured in the horizontal isocenter plane in a diagonal direction with respect to the jaw field.

6.3.4 Measurement of Scatter Factors

The scatter-factors provide information relative to the nominal Linac output. In the Beam Profile Editor the data can be entered arbitrarily normalized.

- The measurement must be performed at the same depth and with the same SSD as for the nominal Linac output measurement (50 mm depth, SSD=1000 mm).
- The scatter factors must be measured for a matrix of combinations of square mMLC fields times square jaw-settings (Figure 3).
- Measure for jaw square field sizes of 6x6, 12x12, 18x18, 42x42, 60x60, 80x80, 100x100 [mm²].
- Measure with mMLC square field sizes of: 6x6, 12x12, 18x18, 24x24, 30x30, 36x36, 42x42, 60x60, 80x80, 100x100 [mm²]
- It is important to measure the fields exactly as indicated.
- It is important to measure all scatter factors where the mMLC size is smaller or equal to the jaw field size (see figure 3, gray fields). If you are sure that the jaw sizes always exceed the mMLC sizes, the scatter factors where mMLC size is larger than jaw size may be set equal to the value where mMLC and jaw size are equal (see figure 3, white fields). In any case, the complete matrix of scatter factors has to be copied to the beam profile editor and no zero values have to remain.

- Note: For the mMLC you can open only 2, 4, 6... leaves, i.e. there is no different choice for the size of the squares (except for adding more square fields to the list above).


mMLC SQUARE FIELD SIZES, mm	JAW SETTINGS, mm						
	6 x 6	12 x 12	18 x 18	42 x 42	60 x 60	80 x 80	100 x 100
6x6	0.605	0.605	0.606	0.607	0.608	0.608	0.609
12x12	0.605	0.756	0.759	0.766	0.766	0.766	0.768
18x18	0.605	0.756	0.796	0.815	0.815	0.815	0.817
24x24	0.605	0.756	0.796	0.840	0.841	0.841	0.842
30x30	0.605	0.756	0.796	0.860	0.860	0.860	0.862
36x36	0.605	0.756	0.796	0.876	0.878	0.880	0.880
42x42	0.605	0.756	0.796	0.888	0.894	0.896	0.897
60x60	0.605	0.756	0.796	0.888	0.931	0.936	0.937
80x80	0.605	0.756	0.796	0.888	0.931	0.969	0.972
100x100	0.605	0.756	0.796	0.888	0.931	0.969	1.000

Figure 3: Example for the table of measured scatter factors. Dark gray fields have to be measured in all cases. White values should be measured, if you plan to use jaw fields which are smaller than mMLC fields. If you do not plan to do so, white fields can be continued using the last measured value in each column, i.e. for a jaw size of 60 x 60 mm² you can take the value of the 60 x 60 mm² mMLC field and copy it downward to all fields for larger mMLC fields.

UNCONTROLLED INFORMATION IF NOT READ ONLINE

6.3.5 Measurement of Diagonal Radial Profiles (Not suited to be entered into BPE directly!)

Radial Factors are dose functions in horizontal direction through the beam axis in several depths. Their purpose is to correct for minor off-axis inhomogeneities. They are only slightly varying functions. Since it is not possible to dismount the mMLC with Novalis they cannot be measured directly. *Radial Profiles* measured in diagonal direction (to reduce boundary effects) contain the same information and can be converted to *Radial Factors*.

The same water phantom setup as for the TMR/PDD measurements can be used. For *Radial Profiles* it is important to have good spatial resolution in the penumbra region (at the diagonal edges). Therefore you have to use a medium size detector with small active volume (pinpoint or smaller). To get results with low noise, it is recommended to move the detector slowly, especially for the lower depth values (to avoid water waves). See figure 4.a. for illustration.

- Open the jaws as well as the mMLC field to 100x100 mm.
- Measure the Radial Profiles in depths of: 5, 14, 25, 50, 100, 200, 350,... [mm]
- Use a radial resolution of at least 2 mm.
- Measure the whole profile ranging from one corner of the field to the other (at least -150 mm to +150 mm in isocenter plane).
- All the profiles have to use the same radius values (either isocentric or real). If necessary use the water phantom software to interpolate.
- Transfer your measurement results to the Excel template.
- Similar to TMR/PDD the Radial Factors/Profiles can be measured in a TMR-like or a PDD-like approach. The radius is always given in the depth where the detector is moving. There is no radius conversion necessary. If radial factors have been measured in a PDD-like approach, the SSD for radial factors must be equal to the SSD used in the PDD measurement. The kind of data measured has to be communicated to BrainLAB HQ.
- Radial Profiles can be normalized arbitrarily.

The fact, that data and the beam profile were sent to BrainLAB does not affect the overall responsibility of the hospital's physicist for the accurateness of the measured values and the correctness of the beam profile!

UNCONTROLLED INFORMATION IF NOT READ ONLINE


Figure 4: Beams eye view indicating the direction of measurement for the radial profiles (a) and the transversal profiles (b). Radial factors (a) are measured in diagonal direction to the jaw field while transversal profiles (b) are measured along and perpendicular to the leaf direction. x-profiles have to be measured with an offset of half the leaf thickness in y-direction (1.5 mm for Novalis) to avoid the inter leaf gap. For y-profiles, either the intraleaf gap or the profile axis has to be shifted by 20 mm from the beam center. In both cases, the detector should be mounted in such a way as to allow maximum spatial resolution along the measuring direction. For Novalis none of these two measurements is suited to be entered into the Beam Profile Editor directly!

UNCONTROLLED INFORMATION IF NOT READ ONLINE

6.4 Additional Measurements

6.4.1 Measurement of dynamic leaf shift (Dynamic IMRT installations, only)

The dynamic leaf shift describes an effective leaf shift due to the round leaf end design of the mMLC. It is determined by measuring the isocenter doses for sliding gaps with different widths (use prepared mMLC files). The value is entered into the beam profile in the "Radiologic Field / Leaf shift dynamic" field.

- Use a large detector (ionization chamber) and position it in the water phantom in a way, that the detector axis is perpendicular to the leaf direction. Set the water surface level to an SSD = 980 mm and adjust the detector at the isocenter (depth = 20 mm).
- Set the jaws to form a square field of 100x100 mm².
- Set the dose rate of the LINAC to 300 MU/min and deliver 300 MU for each field.
- Successively irradiate the dynamic mMLC files "M3_1.d01", "M3_5.d01", ... , "M3_100.d01" (gap sizes: 1, 5, 10, 20, 50, 100 mm) and note the dose values.
- Close the mMLC and measure the leakage dose using the same setting as above (asymmetric gap setting; intra-leaf gap should be 50 mm off the isocenter).
- Set the mMLC to a square field of 100x100 mm² and measure the open field dose using the same settings as above.

 BrainLAB	Working Instruction	WOI 10-26
	Beam Measurement Pencil Beam (Novalis)	Revision 6 page 10 of 10

- The measured dose D can approximately be described by a linear function
 $D - D_{leak} = b(gap + 2\delta) = b \cdot gap + a$, where gap is the nominal gap width (1, 5, ..., 100 mm),
 D_{leak} is the measured mMLC leakage and δ is the effective dynamic leaf shift per leaf. After determination of a and b by linear regression, δ is calculated by $\delta = a/2b$. Please see the Excel sheet "Dynamic Shift" in the Excel template.

6.5 Source Function Correction and Radiological Field

The source function correction is an empirical way of simulating an extended beam source and other effects, that smear out the beam edge. A Gaussian curve with a certain amplitude and a given width (Sigma is a measure for the width) is convolved with the calculated dose distribution. This only influences the penumbra region.

The radiologic field correction allows to correct for small deviations of the radiologic field with respect to the field size defined by the mMLC, which result from gap settings, round leaf end and tongue and groove design. In leaf direction, an offset can be defined by the value "Leaf shift static" in the "Radiologic Field" section of the BPE. Perpendicular to the leaf direction, the offset can be defined by the "Tongue and groove size".

Default values for the Novalis mMLC:

Energy	Source Function Correction		Radiologic Field	
	0 cm depth*	20 cm depth*	Leaf shift static	Tongue and groove size
6 MV	Sigma = 0.18 cm Amplitude = 50%	Sigma = 0.12 cm Amplitude = 50%	0.0 mm	0.0 mm

* Sigma and Amplitude values at all other depths are linearly inter-/extrapolated.

The parameters in the BPE section lightfield correction allow to correct for small deviations of the lightfield. If the lightfield size is correctly adjusted, both values should be set to 0.

Both, the values for the Source Function correction as well as for the Radiologic Field Correction, should not be changed without consulting qualified BrainLAB persona!

7 References

"BrainSCAN", BrainLAB User Manual
 "Beam Profile Editor", BrainLAB User Manual
 WOI Beam Measurement RS and Conformal
 WOI Pencil Beam Measurement (Standard)

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	Working Instruction	WOI 09-10
	Beam Measurements RS	Revision 2 page 1 of 5

1. Purpose

Acquisition of beam data which are suitable for Treatment Planning with BrainSCAN module *Radiosurgery*.

2. Scope

This working instruction is valid for all BrainLAB Support Engineers.

3. Definitions and Abbreviations

OAR : Off-Axis Ratio. Definition: $OAR(c,r,d) = D(c,r,d)/D(c,0,d)$

TMR : Tissue Maximum Ratio. Definition: $TMR(c,d) = D(c,0,d)/D(c,0,d_{max})$

PDD : Percentage Depth-Dose. Converted from TMR by $TMR = \frac{1}{100} * PDD * \left(\frac{SSD + d}{SSD + d_{max}}\right)^2$

S_t : Total Scatter Factor, also referred to as Relative Output Factor.
Definition: $S_t(c) = D(c,r,d_{max})/D(10*10cm^2,0,d_{cat})$

- c : diameter of conical collimator attached to LINAC
- r : off-axis distance perpendicular to central beam axis
- d : depth in water or water-equivalent material
- d_{max} : depth of maximum dose (build-up effect)
- SSD : source-surface distance

UNCONTROLLED INFORMATION IF PRINTED ON PAPER

Originated:	Approved:	Released:
January 26, 1999	January 26, 1999	January 26, 1999
Richard Braun	Stefan Vilsmeier	Per Persson
Product Support Manager	President	Quality Manager


	Working Instruction	WOI 09-10
	Beam Measurements RS	Revision 2 page 2 of 5

4. Instructions

Measurements of TMR, OAR and S_c for all circular stereotactic collimators **MUST** be made with a fixed primary jaw field size setting, usually $5 \times 5 \text{cm}^2$.

OAR measurement set-up

UNCONTROLLED INFORMATION IF PRINTED ON PAPER


OAR can be measured either in a water phantom using ionization chamber or diodes, or with water equivalent material using dosimetry films. When using a water phantom adjust the chamber at isocenter level, 7.5cm below water surface. If you are using an automatic dose measurement system, make a continuous dose measurement in order to get the most accurate and continuous OAR graph. If only a non-automatic system can be used, make sure to measure in 1mm steps in the area of steep dose fall-off, and every 3-5 mm at the flat parts of the dose curve until the dose drops to less than 2% of the maximum dose.


Make sure to measure the profile in positive and negative range. According to the definition of OAR, normalize the data to the dose measured on central beam axis $r=0$, and average corresponding data points from the negative and positive branch.

OAR has to be measured for each collimator separately using the same jaw position.

TMR measurement set-up

PDD measurement set-up

UNCONTROLLED INFORMATION IF PRINTED ON PAPER

Usually, this is done with ionization chambers only. The chamber is fixed inside the water phantom at isocenter position. At the beginning, the device is placed just on the water surface. Then water is filled up to a level of 25cm above isocenter. During filling, the LINAC is continuously irradiating and the chamber is measuring the dose. In case the measurement system can not pick up the data continuously, BrainLAB recommends to use the following water heights for reconstructing the TMR-profile (for 6MV beam): 0, 5, 8, 12, 14, 15, 16, 18, 20, 25, 30, 40, 50, 75, 100, 125, 150, 175, 200, 225, 250mm. For photon beams with higher energy the depths should be modified corresponding to the expected changes in depth-dose curve.

Adjust water surface to d_{max} above isocenter ($SSD = 100 - d_{max}$). Position chamber at water surface. Then move the chamber 25cm down into the water. Keep the LINAC continuously irradiating and the chamber measuring during movement. In case the measurement system can not pick up the data continuously, BrainLAB recommends to use the following depths for reconstructing the PDD-profile (for 6MV beam e.g.): 0, 5, 8, 12, 14, 15, 16, 18, 20, 25, 30, 40, 50, 75, 100, 125, 150, 175, 200, 225, 250mm. For photon beams with higher energy the depths should be modified corresponding to the expected changes in depth-dose curve.

Depth-doses show a maximum dose at a particular depth that depends, beneath others, on collimator size. But the algorithms used for dosage calculation require all data to be normalized to the same depth. This depth will be referred to as d_{max} . Following this normalization instruction, it might happen that some graphs exceed unity.

TMR is defined as:

$$TMR(c, d) = D(c, 0, d) / D(c, 0, d_{max})$$

- c : diameter of cone at isocenter
- r : off-axis distance
- d : depth in tissue


To get data comparable to TMR the data measured with the PDD method need to be converted to TMR using:

$$TMR(c, d) = PDD(c, d, SSD) \cdot \left(\frac{SSD + d}{SSD + d_{max}} \right)^2$$


The depth-dose data has to be measured for each collimator separately using the same jaw position.

Scatter Factor S_t (Relative Output Factor)

$S_c(c, 0, d_{max})$ measurement set-up

$S_c(c, 0, d_{cal})$ measurement set-up
(Hospital calibration setup, e.g.)

The scatter factor accounts for the attenuation of the beams that is due to the focusing effect of the collimator in comparison to an "open field" irradiation.

Measure the dose at the reference point d_{max} ($SSD = 100 - d_{max}$), i.e. the depth where the PDD/TMR have been normalized. Apply the same dose that is delivered when measuring the nominal LINAC output and when calibrating the LINAC, e.g. $MU_{cal} = 100$ MU. Make three measurements (for each collimator) and take the average to account for inaccurate dose delivery/measurement. Keep the jaws in the same position as for the measurements of PDD/TMR and OAR.

	Working Instruction	WOI 09-10
	Beam Measurements RS	Revision 2 page 5 of 5

Then measure the dose for the hospital calibration setup (e.g. SID=SSD=100cm, $d_{max}=d_{cal}$, fieldsize=10x10cm²) at depth d_{cal} . Again, take the average out of three measurements.

S_t is defined as $S_{tc} = D(c, 0, d_{max}) / D(10 \times 10 \text{cm}^2, 0, d_{cal})$

Reason: Measuring at hospital calibration setup gives an electrometer reading of M_{cal} that is according to the hospital a dose D_{cal} (e.g. 1Gy). Thus we have calibrated the ionization chamber and calculate a dose from given electrometer readouts.

S_t has to be measured for each collimator separately using the same jaw position.

5. References

"Beam Shape Editor", BrainLAB Users Manual

"Dose Algorithm", BrainLAB Users Manual

Rice R K, Hansen J L, Svensson G K, Siddon R L; *Measurements of dose distributions in small 6MV x-rays.* (1987) Phys. Med. Biol. 32 pp1087-1099.

UNCONTROLLED INFORMATION IF PRINTED ON PAPER

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

Various build-up caps for compact chambers and farmer type chambers in PMMA or other materials on request.
Farmer type chambers: For Cobalt energies the protection cap can be used as build-up caps.

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Specifications	Item 300-672
PTW	31014
Measuring Volume:	0.015 cm ³ Universal
Nominal Useful Range:	⁶⁰ Co to 50 MV 6 to 50 MeV
Photons:	
Electrons:	
Response:	400 pC/Gy
Leakage Current:	<E ± 4 fA
Polarizing Voltage Effect:	< 1%
Cable Leakage:	± 1 pC/(Gy • cm)
Wall Material & Thickness:	0.57 mm PMMA 0.09 mm Graphite
Wall Density:	PMMA: 1.19 g/cm ³ Graphite: 1.84 g/cm ³
Area Density:	84.4 mg/cm ²
Electrode Material:	
Ion Collection Time:	At 99.98 R 400 V: 20 μs